

产品名称: KIR2DL2/3 Rabbit Polyclonal Antibody
产品货号: APRab00620



产品概述 (Summary)

产品名称 (Production Name)	KIR2DL2/3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KIR2DL3/KIR2DS2
别名 (Alternative Names)	KIR2DL3; CD158B2; KIRCL23; NKAT2; Killer cell immunoglobulin-like receptor 2DL3; CD158 antigen-like family member B2; KIR-023GB; Killer inhibitory receptor cl 2-3; MHC class I NK cell receptor; NKAT2a; NKAT2bNatural killer-associated transcript 2; NKAT-2; p58 natural killer cell receptor clone CL-6; p58 NK receptor CL-6; p58.2 MHC class-I-specific NK receptor; CD158b2; KIR2DS2; CD158J; NKAT5; Killer cell immunoglobulin-like receptor 2DS2
基因 ID (Gene ID)	3804/100132285
蛋白 ID (SwissProt ID)	P43628/P43631.

产品应用 (Application)

产品名称: KIR2DL2/3 Rabbit Polyclonal Antibody
产品货号: AP Rab00620



稀释比 (Dilution Ratio)

WB 1:500-1:1000, ELISA 1:5000-1:20000

蛋白分子量 (Molecular Weight)

Calculated MW: 38 kDa; Observed MW: 38 kDa


研究背景 (Background)

Receptor on natural killer (NK) cells for HLA-C alleles (HLA-Cw1, HLA-Cw3 and HLA-Cw7). Inhibits the activity of NK cells thus preventing cell lysis.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of KIR2DL2/3 in HeLa lysates using KIR2DL2/3 antibody.

注意事项 (Note)

For research use only .