

产品名称: Lymphotactin Rabbit Polyclonal Antibody
产品货号: APRab00621



产品概述 (Summary)

产品名称 (Production Name)	Lymphotactin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	XCL1/XCL2
别名 (Alternative Names)	XCL1; LTN; SCYC1; Lymphotactin; ATAC; C motif chemokine 1; Cytokine SCM-1; Lymphotaxin; SCM-1-alpha; Small-inducible cytokine C1; XC chemokine ligand 1; XCL2; SCYC2; Cytokine SCM-1 beta; C motif chemokine 2; XC chemokine ligand 2
基因 ID (Gene ID)	6375/6846
蛋白 ID (SwissProt ID)	P47992/Q9UBD3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 13 kDa; Observed MW: 13 kDa

产品名称: Lymphotactin Rabbit Polyclonal Antibody
产品货号: APRab00621




研究背景 (Background)

Chemotactic activity for lymphocytes but not for monocytes or neutrophils. In thymus, mediates medullary accumulation of thymic dendritic cells and contributes to regulatory T cell development, playing a role in self-tolerance establishment.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of Lymphotactin in K562 lysates using Lymphotactin antibody.

注意事项 (Note)

For research use only .