

产品名称: CD85f Rabbit Polyclonal Antibody
产品货号: APRab00623



产品概述 (Summary)

产品名称 (Production Name)	CD85f Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LILRA5
别名 (Alternative Names)	LILRA5; ILT11; LILRB7; LIR9; Leukocyte immunoglobulin-like receptor subfamily A member 5; CD85 antigen-like family member F; Immunoglobulin-like transcript 11; ILT-11; Leukocyte immunoglobulin-like receptor 9; LIR-9; CD85f
基因 ID (Gene ID)	353514
蛋白 ID (SwissProt ID)	A6NI73.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 33 kDa

产品名称: CD85f Rabbit Polyclonal Antibody
产品货号: APRab00623




研究背景 (Background)

May play a role in triggering innate immune responses.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD85f in HeLa lysates using CD85f antibody.

注意事项 (Note)

For research use only .