

产品名称: DCL-1 Rabbit Polyclonal Antibody
产品货号: APRab00624



产品概述 (Summary)

产品名称 (Production Name)	DCL-1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD302 DCL-1; CD302; CLEC13A; DCL1; KIAA0022; CD302 antigen; C-type lectin
别名 (Alternative Names)	BIMLEC; C-type lectin domain family 13 member A; DEC205-associated C-type lectin 1; Type I transmembrane C-type lectin receptor DCL-1; CD302
基因 ID (Gene ID)	9936
蛋白 ID (SwissProt ID)	Q8IX05.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 26 kDa

产品名称: DCL-1 Rabbit Polyclonal Antibody
产品货号: APRab00624




研究背景 (Background)

Potential multifunctional C-type lectin receptor that may play roles in endocytosis and phagocytosis as well as in cell adhesion and migration.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of DCL-1 in HeLa lysates using DCL-1 antibody.

注意事项 (Note)

For research use only .