

产品名称: Elongation Factor 1A1/2/P5 Rabbit Polyclonal Antibody
产品货号: APRab00625



产品概述 (Summary)

产品名称 (Production Name)	Elongation Factor 1A1/2/P5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EEF1A1/EEF1A2/EEF1A1P5 EEF1A1; EEF1A; EF1A; LENG7; Elongation factor 1-alpha 1; EF-1-alpha-1; Elongation factor Tu; EF-Tu; Eukaryotic elongation factor 1 A-1; eEF1A-1; Leukocyte receptor cluster member 7; EEF1A2; EEF1AL; STN; Elongation factor 1-alpha 2; EF-1-alpha-2; Eukaryotic elongation factor 1 A-2; eEF1A-2; Statin-S1; Putative elongation factor 1-alpha-like 3; EF-1-alpha-like 3; Eukaryotic elongation factor 1 A-like 3; eEF1A-like 3
别名 (Alternative Names)	
基因 ID (Gene ID)	1915/1917/158078
蛋白 ID (SwissProt ID)	P68104/Q05639/Q5VTE0.

产品应用 (Application)

产品名称: Elongation Factor 1A1/2/P5 Rabbit Polyclonal Antibody
产品货号: APRab00625



稀释比 (Dilution Ratio)

WB 1:500-1:1000, ELISA 1:5000-1:20000

蛋白分子量 (Molecular Weight)

Calculated MW: 50 kDa; Observed MW: 50 kDa


研究背景 (Background)

This protein promotes the GTP-dependent binding of aminoacyl-tRNA to the A-site of ribosomes during protein biosynthesis.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Elongation Factor 1A1/2/P5 in HepG2 lysates using Elongation Factor 1A1/2/P5 antibody.

注意事项 (Note)

For research use only .