

产品名称: FGF11 Rabbit Polyclonal Antibody
产品货号: APRab00628



产品概述 (Summary)

产品名称 (Production Name)	FGF11 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGF11
别名 (Alternative Names)	FGF11; FHF3; Fibroblast growth factor 11; FGF-11; Fibroblast growth factor homologous factor 3; FHF-3
基因 ID (Gene ID)	2256
蛋白 ID (SwissProt ID)	Q92914.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 25 kDa

研究背景 (Background)

产品名称: FGF11 Rabbit Polyclonal Antibody
产品货号: APRab00628




Probably involved in nervous system development and function.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of FGF11 in NIH3T3, HeLa lysates using FGF11 antibody.

注意事项 (Note)

For research use only .