

产品名称: AChR alpha1 Rabbit Polyclonal Antibody
产品货号: APRab00629



产品概述 (Summary)

产品名称 (Production Name)	AChR alpha1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHRNA1
别名 (Alternative Names)	CHRNA1; ACHRA; CHNRA; Acetylcholine receptor subunit alpha
基因 ID (Gene ID)	1134
蛋白 ID (SwissProt ID)	P02708.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 55 kDa

研究背景 (Background)

Upon acetylcholine binding, the AChR responds by an extensive change in conformation that affects all subunits and leads

产品名称: AChR alpha1 Rabbit Polyclonal Antibody
产品货号: AP Rab00629




to opening of an ion-conducting channel across the plasma membrane.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of AChR alpha1 in NIH3T3 lysates using AChR alpha1 antibody.

注意事项 (Note)

For research use only .