

产品名称: STRA8 Rabbit Polyclonal Antibody
产品货号: APRab00630



产品概述 (Summary)

产品名称 (Production Name)	STRA8 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STRA8
别名 (Alternative Names)	STRA8; Stimulated by retinoic acid gene 8 protein homolog
基因 ID (Gene ID)	346673
蛋白 ID (SwissProt ID)	Q7Z7C7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 66 kDa

研究背景 (Background)

This gene encodes a retinoic acid-responsive protein. A homologous protein in mouse has been shown to be involved in

产品名称: STRA8 Rabbit Polyclonal Antibody
产品货号: APRab00630




the regulation of meiotic initiation in both spermatogenesis and oogenesis, though feature differences between the mouse and human proteins suggest that these homologs are not entirely functionally equivalent. It is thought that this gene may play a role in spermatogenesis in humans.

研究领域 (Research Area)

Stem Cells

图片 (Image Data)



Western blot analysis of STRA8 in NIH3T3, 293 lysates using STRA8 antibody.

注意事项 (Note)

For research use only .