

产品名称: PI3 Kinase C2 gamma Rabbit Polyclonal Antibody
产品货号: APRab00639



产品概述 (Summary)

产品名称 (Production Name)	PI3 Kinase C2 gamma Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PIK3C2G
别名 (Alternative Names)	PIK3C2G; Phosphatidylinositol 4-phosphate 3-kinase C2 domain-containing subunit gamma; PI3K-C2-gamma; PtdIns-3-kinase C2 subunit gamma; Phosphoinositide 3-kinase-C2-gamma
基因 ID (Gene ID)	5288
蛋白 ID (SwissProt ID)	O75747.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 166 kDa; Observed MW: 160 kDa

产品名称: PI3 Kinase C2 gamma Rabbit Polyclonal Antibody
产品货号: APRab00639




研究背景 (Background)

Generates phosphatidylinositol 3-phosphate (PtdIns3P) and phosphatidylinositol 3,4-bisphosphate (PtdIns(3,4)P₂) that act as second messengers (By similarity).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of PI3 Kinase C2 gamma in A549 lysates using PI3 Kinase C2 gamma antibody.

注意事项 (Note)

For research use only .