

**产品名称: Cripto 3 Rabbit Polyclonal Antibody**  
**产品货号: APRab00641**



## 产品概述 (Summary)

|                        |                                     |
|------------------------|-------------------------------------|
| 产品名称 (Production Name) | Cripto 3 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody          |
| 宿主 (Host)              | Rabbit                              |
| 应用 (Application)       | WB,ELISA                            |
| 种属反应性 (Reactivity)     | Human                               |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Unmodified                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                            |
| 克隆 (Clonality)      | Polyclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Affinity Purification                                                                          |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                                                                            |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | TDGF1P3<br>TDGF1P3; CRIPTO3; TDGF2; TDGF3; Putative teratocarcinoma-derived growth factor 3; Cripto-3 growth factor; Epidermal growth factor-like cripto protein CR3; Teratocarcinoma-derived growth factor 1 pseudogene 3 |
| 别名 (Alternative Names) |                                                                                                                                                                                                                            |
| 基因 ID (Gene ID)        | 6998.0                                                                                                                                                                                                                     |
| 蛋白 ID (SwissProt ID)   | P51864.                                                                                                                                                                                                                    |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,ELISA 1:5000-1:20000       |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 21 kDa; Observed MW: 21 kDa |

**产品名称: Cripito 3 Rabbit Polyclonal Antibody**  
**产品货号: APRab00641**




## 研究背景 (Background)

Could play a role in the determination of the epiblastic cells that subsequently give rise to the mesoderm.

## 研究领域 (Research Area)

Cell Biology

## 图片 (Image Data)



Western blot analysis of Cripito 3 in HepG2, 293T, AD293, HeLa lysates using Cripito 3 antibody.

## 注意事项 (Note)

For research use only.