

产品名称: Collagen XII alpha 1 Rabbit Polyclonal Antibody
产品货号: APRab00650



产品概述 (Summary)

产品名称 (Production Name)	Collagen XII alpha 1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	COL12A1
别名 (Alternative Names)	Collagen alpha-1(XII) chain
基因 ID (Gene ID)	1303
蛋白 ID (SwissProt ID)	Q99715.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

Type XII collagen interacts with type I collagen-containing fibrils, the COL1 domain could be associated with the surface of

产品名称: Collagen XII alpha 1 Rabbit Polyclonal Antibody
产品货号: APRab00650




the fibrils, and the COL2 and NC3 domains may be localized in the perifibrillar matrix.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human muscle using Collagen XII alpha 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .