

产品名称: Cytokeratin Rabbit Polyclonal Antibody
产品货号: AP Rab00651



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Keratin KRT2A; KRT2E; KRT2B; KRT2P; K6A; KRT6D; K6B; KRTL1; KRT6E; K6IRS1; KB34;
别名 (Alternative Names)	KRT6IRS1; K6IRS2; KB35; KRT6; KRT6IRS2; K6IRS3; KB36; KRT6IRS3; K6IRS4; KB37; KRT5C; KRT6IRS4; K6HF; KB18; K6L; KB38; KRT6L; SCL; CYK8; KRTHB4
基因 ID (Gene ID)	3849.0
蛋白 ID (SwissProt ID)	P35908/Q01546/P12035/P13647/P02538/P04259/P48668/Q3SY84/ Q14CN4/Q86Y46/Q7RTS7/O95678/Q5XKE5/P08729/P05787/Q9NSB2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 65 kDa; Observed MW: 70 kDa

产品名称: Cytokeratin Rabbit Polyclonal Antibody
产品货号: APRab00651




研究背景 (Background)

Probably contributes to terminal cornification. Associated with keratinocyte activation, proliferation and keratinization.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Cytokeratin in mouse lung, mouse brain, mouse kidney, mouse heart lysates using Cytokeratin antibody.

注意事项 (Note)

For research use only .