

产品名称: TNF alpha Rabbit Polyclonal Antibody
产品货号: APRab00662



产品概述 (Summary)

产品名称 (Production Name)	TNF alpha Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNF
别名 (Alternative Names)	TNF; TNFA; TNFSF2; Tumor necrosis factor; Cachectin; TNF-alpha; Tumor necrosis factor ligand superfamily member 2; TNF-a
基因 ID (Gene ID)	7124
蛋白 ID (SwissProt ID)	P01375.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 26 kDa

研究背景 (Background)

产品名称: TNF alpha Rabbit Polyclonal Antibody
产品货号: APRab00662




This gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR2. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, and cancer. Knockout studies in mice also suggested the neuroprotective function of this cytokine.


研究领域 (Research Area)

Immunology

图片 (Image Data)




Western blot analysis of TNF alpha in mouse liver lysates using TNF alpha antibody.




Immunohistochemistry analysis of paraffin-embedded Human heart using TNF alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.


产品名称: TNF alpha Rabbit Polyclonal Antibody
产品货号: APRab00662



Immunohistochemistry analysis of paraffin-embedded Human heart using TNF alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemistry analysis of paraffin-embedded Human colon using TNF alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human colon using TNF alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: TNF alpha Rabbit Polyclonal Antibody
产品货号: APRab00662



Immunohistochemistry analysis of paraffin-embedded Human heart using TNF alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .