

**产品名称: Phospho-Src (Tyr419) Rabbit Polyclonal Antibody**  
**产品货号: AP Rab00684**



## 产品概述 (Summary)

|                        |                                                 |
|------------------------|-------------------------------------------------|
| 产品名称 (Production Name) | Phospho-Src (Tyr419) Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody                      |
| 宿主 (Host)              | Rabbit                                          |
| 应用 (Application)       | WB,FC                                           |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                                 |

## 产品性能 (Performance)

|                     |                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                       |
| 修饰 (Modification)   | Phosphorylated                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                                |
| 克隆 (Clonality)      | Polyclonal                                                                                         |
| 形式 (Form)           | Liquid                                                                                             |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.           |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Chromatography                                                                            |

## 免疫原信息 (Immunogen)

|                        |                                                                                                 |
|------------------------|-------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | SRC                                                                                             |
| 别名 (Alternative Names) | SRC; SRC1; Proto-oncogene tyrosine-protein kinase Src; Proto-oncogene c-Src; pp60c-src; p60-Src |
| 基因 ID (Gene ID)        | 6714                                                                                            |
| 蛋白 ID (SwissProt ID)   | P12931.                                                                                         |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,FC 1:50-1:100              |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 60 kDa; Observed MW: 60 kDa |

## 研究背景 (Background)

**产品名称: Phospho-Src (Tyr419) Rabbit Polyclonal Antibody**  
**产品货号: AP Rab00684**




Src activity is regulated by tyrosine phosphorylation at two sites, but with opposing effects. While phosphorylation at Tyr416 in the activation loop of the kinase domain upregulates enzyme activity, phosphorylation at Tyr527 in the carboxy-terminal tail by Csk renders the enzyme less active (2).

## 研究领域 (Research Area)

Signal Transduction

## 图片 (Image Data)



Western blot analysis of Phospho-Src (Y419) in A431 lysates treated with pervanadate using Phospho-Src (Tyr419) antibody.

## 注意事项 (Note)

For research use only .