

**产品名称:** Phospho-IKB alpha (Ser32) Rabbit Polyclonal Antibody

**产品货号:** APRab00694

### 产品概述 (Summary)

|                   |                                                                                                    |
|-------------------|----------------------------------------------------------------------------------------------------|
| 描述(Description)   | Rabbit polyclonal Antibody                                                                         |
| 宿主(Host)          | Rabbit                                                                                             |
| 应用(Application)   | WB,IP                                                                                              |
| 种属反应(Reactivity)  | Human                                                                                              |
| 偶联物(Conjugation)  | Unconjugated                                                                                       |
| 修饰(Modification)  | Phosphorylated                                                                                     |
| 同种型(Isotype)      | IgG                                                                                                |
| 克隆(Clonality)     | Polyclonal                                                                                         |
| 剂型(Form)          | Liquid                                                                                             |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.           |
| 储存溶液(Buffer)      | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| 纯化 (Purification) | Affinity Chromatography                                                                            |

### 产品应用(Application)

|                       |                                            |
|-----------------------|--------------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:500-1:1000,IP 1:20-1:50               |
| 分子量(Molecular Weight) | Calculated MW: 36 kDa; Observed MW: 39 kDa |

### 抗原信息(Antigen Information)

|                       |                                                                                                                                                                 |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名(Gene Name)        | NFKBIA                                                                                                                                                          |
| 别名(Alternative Names) | NFKBIA; IKBA; MAD3; NFKBI; NF-kappa-B inhibitor alpha; I-kappa-B-alpha; Ikb-alpha; IkappaBalpha; Major histocompatibility complex enhancer-binding protein MAD3 |
| 基因 ID(Gene ID)        | 4792                                                                                                                                                            |
| SwissProt ID          | P25963                                                                                                                                                          |
| 免疫原(Immunogen)        | A synthetic Phosphorylated peptide corresponding to residues target protein                                                                                     |

### 研究背景 (Background)

NFKB1 (MIM 164011) or NFKB2 (MIM 164012) is bound to REL (MIM 164910), RELA (MIM 164014), or RELB (MIM 604758) to form the NFKB complex. The NFKB complex is inhibited by I-kappa-B proteins (NFKBIA or NFKBIB, MIM 604495), which

inactivate NF-kappa-B by trapping it in the cytoplasm.


### 研究领域 (Research Area)

Epigenetics and Nuclear Signaling

### 注意事项 (Note)

For research use only.

### 图片 (Image Data)



Western blot analysis of Phospho-IKB alpha (S32) in HeLa lysates treated with Calyculin A and TNFa using Phospho-IKB alpha (Ser32) antibody.

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>