

产品名称: Phospho-CDC37 (Ser13) Rabbit Polyclonal Antibody
产品货号: APRab00707



产品概述 (Summary)

产品名称 (Production Name)	Phospho-CDC37 (Ser13) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDC37
别名 (Alternative Names)	CDC37; CDC37A; Hsp90 co-chaperone Cdc37; Hsp90 chaperone protein kinase-targeting subunit; p50Cdc37
基因 ID (Gene ID)	11140
蛋白 ID (SwissProt ID)	Q16543.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 50 kDa

产品名称: Phospho-CDC37 (Ser13) Rabbit Polyclonal Antibody
产品货号: APRab00707




研究背景 (Background)

CDC37 is an important component of the HSP90 chaperone complex. It was initially identified for its involvement in cell-cycle progression and was later found to have a much broader role as a chaperone for a wide variety of kinases and other proteins. CDC37 protein has an amino-terminal kinase binding domain followed by a central HSP90 binding domain.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Phospho-CDC37 (S13) in (1) Jurkat lysates treated with Alkaline Phosphatase; (2) Jurkat lysates using Phospho-CDC37 (Ser13) antibody.

注意事项 (Note)

For research use only .