

产品名称: Phospho-BRAF (Thr401) Rabbit Polyclonal Antibody
产品货号: APRab00709



产品概述 (Summary)

产品名称 (Production Name)	Phospho-BRAF (Thr401) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	BRAF
别名 (Alternative Names)	BRAF; BRAF1; RAFB1; Serine/threonine-protein kinase B-raf; Proto-oncogene B-Raf; p94; v-Raf murine sarcoma viral oncogene homolog B1
基因 ID (Gene ID)	673
蛋白 ID (SwissProt ID)	P15056.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 84 kDa; Observed MW: 84 kDa

产品名称: Phospho-BRAF (Thr401) Rabbit Polyclonal Antibody
产品货号: APRab00709




研究背景 (Background)

Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron. Frequently mutated in thyroid cancers, skin melanomas and at lower frequency in a wide range of human cancers. An activating mutation, mimicking phosphorylation of the activation loop, is seen in 60% of malignant melanoma samples. Raf mutations are generally exclusive to Ras activating mutations.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Phospho-B Raf (T401) in (1) PC-12 lysates; (2) PC-12 lysates treated with TPA using Phospho-BRAF (Thr401) antibody.

注意事项 (Note)

For research use only .