

产品名称: Phospho-SIRT1 (Ser47) Rabbit Polyclonal Antibody
产品货号: APRab00715



产品概述 (Summary)

产品名称 (Production Name)	Phospho-SIRT1 (Ser47) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SIRT1
别名 (Alternative Names)	SIRT1; SIR2L1; NAD-dependent protein deacetylase sirtuin-1; hSIRT1; Regulatory protein SIR2 homolog 1; SIR2-like protein 1; hSIR2
基因 ID (Gene ID)	23411
蛋白 ID (SwissProt ID)	Q96EB6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 82 kDa; Observed MW: 120 kDa

产品名称: Phospho-SIRT1 (Ser47) Rabbit Polyclonal Antibody
产品货号: APRab00715




研究背景 (Background)

This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Phospho-SIRT1(Ser47) in A549 lysates using Phospho-SIRT1 (Ser47) antibody.

注意事项 (Note)

For research use only .