

**产品名称: Phospho-Neurofilament Medium (Ser614/Ser619)
Rabbit Polyclonal Antibody**
产品货号: APRab00716



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Neurofilament Medium (Ser614/Ser619) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEFM
别名 (Alternative Names)	NEFM; NEF3; NFM; Neurofilament medium polypeptide; NF-M; 160 kDa neurofilament protein; Neurofilament 3; Neurofilament triplet M protein
基因 ID (Gene ID)	4741
蛋白 ID (SwissProt ID)	P07197.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 102 kDa; Observed MW: 160 kDa

产品名称: Phospho-Neurofilament Medium (Ser614/Ser619)
Rabbit Polyclonal Antibody
产品货号: APRab00716


研究背景 (Background)

Integrins are heterodimers composed of noncovalently associated transmembrane α and β subunits. The 16 α and 8 β subunits heterodimerize to produce more than 20 different receptors. Antibodies to this protein are useful to identify neurons and their processes in tissue sections and in tissue culture.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Phospho-Neurofilament Medium (Ser614/619) in K562, Ramos, Daudi lysates using Phospho-Neurofilament Medium (Ser614/Ser619) antibody.

注意事项 (Note)

For research use only .