

产品名称: Phospho-MEK2 (Thr394) Rabbit Polyclonal Antibody
产品货号: APRab00791



产品概述 (Summary)

产品名称 (Production Name)	Phospho-MEK2 (Thr394) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAP2K2
别名 (Alternative Names)	MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity mitogen-activated protein kinase kinase 2; MAP kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2; MEK 2
基因 ID (Gene ID)	5605
蛋白 ID (SwissProt ID)	P36507.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 44 kDa

产品名称: Phospho-MEK2 (Thr394) Rabbit Polyclonal Antibody
产品货号: APRab00791




研究背景 (Background)

Catalyzes the concomitant phosphorylation of a threonine and a tyrosine residue in a Thr-Glu-Tyr sequence located in MAP kinases. Activates the ERK1 and ERK2 MAP kinases.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Phospho-MEK2 (Thr394) in various lysates using Phospho-MEK2 (Thr394) antibody.



Western blot analysis of Phospho-MEK2 (Thr394) in VEC lysates using Phospho-MEK2 (T394) antibody.

注意事项 (Note)

For research use only .