

产品名称: Phospho-CaMKII (Thr305) Rabbit Polyclonal Antibody
产品货号: APRab00792



产品概述 (Summary)

产品名称 (Production Name)	Phospho-CaMKII (Thr305) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAMK2A
别名 (Alternative Names)	CAMK2A; CAMKA; KIAA0968; Calcium/calmodulin-dependent protein kinase type II subunit alpha; CaM kinase II subunit alpha; CaMK-II subunit alpha; CAMK2B; CAM2; CAMK2; CAMKB; Calcium/calmodulin-dependent protein kinase type II subunit beta; Ca
基因 ID (Gene ID)	815/816/817
蛋白 ID (SwissProt ID)	Q9UQM7/Q13554/Q13557.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
----------------------	---

产品名称: Phospho-CaMKII (Thr305) Rabbit Polyclonal Antibody
产品货号: APRab00792



蛋白分子量 (Molecular Weight) Calculated MW: 54 kDa; Observed MW: 54 kDa


研究背景 (Background)

CaM-kinase II (CAMK2) is a prominent kinase in the central nervous system that may function in long-term potentiation and neurotransmitter release. Member of the NMDAR signaling complex in excitatory synapses it may regulate NMDAR-dependent potentiation of the AMPAR and synaptic plasticity.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Phospho-CaMKII (Thr305) in various lysates using Phospho-CaMKII (Thr305) antibody.



Western blot analysis of Phospho-CaMKII (Thr305) in 453 lysates using Phospho-CaMKIIα/β/δ (T305) antibody.

注意事项 (Note)

For research use only .