

产品名称: LEFTY Rabbit Polyclonal Antibody
产品货号: APRab00796



产品概述 (Summary)

产品名称 (Production Name)	LEFTY Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LEFTY1/LEFTY2 LEFTY1; LEFTB; LEFTYB; Left-right determination factor 1; Left-right determination factor B; Protein lefty-1; Protein lefty-B; LEFTY2; EBAF; LEFTA;
别名 (Alternative Names)	LEFTYA; TGFB4; PSEC0024; Left-right determination factor 2; Endometrial bleeding-associated factor; Left-right determination factor A; Protein lefty-2; Protein lefty-A; Transforming growth factor beta-4; TGF-beta-4
基因 ID (Gene ID)	10637/7044
蛋白 ID (SwissProt ID)	O75610/O00292.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
----------------------	--------------------------------------

产品名称: LEFTY Rabbit Polyclonal Antibody
产品货号: APRab00796



蛋白分子量 (Molecular Weight) Calculated MW: 41 kDa; Observed MW: 41 kDa


研究背景 (Background)

Required for left-right axis determination as a regulator of LEFTY2 and NODAL.


研究领域 (Research Area)

Stem Cells

图片 (Image Data)



Western blot analysis of LEFTY in HepG2UV lysates using LEFTY antibody.



Western blot analysis of LEFTY in HEPG2UV lysates using Lefty antibody.

注意事项 (Note)

For research use only .