

产品名称: TNFAIP8L3 Rabbit Polyclonal Antibody
产品货号: APRab00801



产品概述 (Summary)

产品名称 (Production Name)	TNFAIP8L3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNFAIP8L3
别名 (Alternative Names)	Tumor necrosis factor alpha-induced protein 8-like protein 3; TNF alpha-induced protein 8-like protein 3; TNFAIP8-like protein 3
基因 ID (Gene ID)	388121
蛋白 ID (SwissProt ID)	Q5GJ75.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 33 kDa

研究背景 (Background)

产品名称: TNFAIP8L3 Rabbit Polyclonal Antibody
产品货号: APRab00801




Acts as a lipid transfer protein.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of TNFAIP8L3 in 293T , mouse brain lysates using TNFAIP8L3 antibody.

注意事项 (Note)

For research use only .