

产品名称: Integrin beta 7 Rabbit Polyclonal Antibody
产品货号: APRab00802



产品概述 (Summary)

产品名称 (Production Name)	Integrin beta 7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITGB7
别名 (Alternative Names)	Integrin beta-7 (Gut homing receptor beta subunit)
基因 ID (Gene ID)	3695
蛋白 ID (SwissProt ID)	P26010.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 87 kDa; Observed MW: 130 kDa

研究背景 (Background)

Integrin alpha-4/beta-7 (Peyer patches-specific homing receptor LPAM-1) is an adhesion molecule that mediates

产品名称: Integrin beta 7 Rabbit Polyclonal Antibody
产品货号: AP Rab00802




lymphocyte migration and homing to gut-associated lymphoid tissue (GALT).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Integrin beta 7 in mouse heart lysates using Integrin beta 7 antibody.

注意事项 (Note)

For research use only .