

产品名称: MAP4K2/3 Rabbit Polyclonal Antibody
产品货号: APRab00806



产品概述 (Summary)

产品名称 (Production Name)	MAP4K2/3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAP4K2/MAP4K3 Mitogen-activated protein kinase kinase kinase kinase 2; B lymphocyte serine/threonine-protein kinase; Germinal center kinase; GC kinase;
别名 (Alternative Names)	MAPK/ERK kinase kinase kinase 2; MEK kinase kinase 2; MEKKK 2; Rab8-interacting protein; Mitogen-activated protein kinase kinase kinase kinase 3; Germinal center kinase-related protein kinase; GLK; MAPK/ERK kinase kinase kinase 3; MEK kinase kinase 3; MEKKK 3
基因 ID (Gene ID)	5871/8491
蛋白 ID (SwissProt ID)	Q12851/Q8IVH8.

产品应用 (Application)

产品名称: MAP4K2/3 Rabbit Polyclonal Antibody
产品货号: APRab00806



稀释比 (Dilution Ratio)

WB 1:500-1:1000, ELISA 1:5000-1:20000

蛋白分子量 (Molecular Weight)

Calculated MW: 92 kDa; Observed MW: 100 kDa


研究背景 (Background)

Serine/threonine-protein kinase which acts as an essential component of the MAP kinase signal transduction pathway. Acts as a MAPK kinase kinase (MAP4K) and is an upstream activator of the stress-activated protein kinase/c-Jun N-terminal kinase (SAP/JNK) signaling pathway and to a lesser extent of the p38 MAPKs signaling pathway.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of MAP4K2/3 in 293T, hela lysates using MAP4K2/3 antibody.

注意事项 (Note)

For research use only .