

**产品名称: Phospho-JAK1 (Tyr1022) Rabbit Polyclonal Antibody**  
**产品货号: APRab00820**



## 产品概述 (Summary)

|                        |                                                   |
|------------------------|---------------------------------------------------|
| 产品名称 (Production Name) | Phospho-JAK1 (Tyr1022) Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody                        |
| 宿主 (Host)              | Rabbit                                            |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA                               |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                                   |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Phosphorylated                                                                                 |
| 同种型 (Isotype)       | IgG                                                                                            |
| 克隆 (Clonality)      | Polyclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Affinity Purification                                                                          |

## 免疫原信息 (Immunogen)

|                        |                                                                         |
|------------------------|-------------------------------------------------------------------------|
| 基因名 (Gene Name)        | JAK1                                                                    |
| 别名 (Alternative Names) | JAK1; JAK1A; JAK1B; Tyrosine-protein kinase JAK1; Janus kinase 1; JAK-1 |
| 基因 ID (Gene ID)        | 3716                                                                    |
| 蛋白 ID (SwissProt ID)   | P23458.                                                                 |

## 产品应用 (Application)

|                          |                                                                       |
|--------------------------|-----------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 133 kDa; Observed MW: 133 kDa                          |

## 研究背景 (Background)

**产品名称: Phospho-JAK1 (Tyr1022) Rabbit Polyclonal Antibody**  
**产品货号: APRab00820**




Tyrosine kinase of the non-receptor type, involved in the IFN-alpha/beta/gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

## 研究领域 (Research Area)

Cell Biology

## 图片 (Image Data)



Western blot analysis of Phospho-JAK1 (Tyr1022) in various lysates using Phospho-JAK1 (Tyr1022) antibody.

## 注意事项 (Note)

For research use only .