

产品名称: Phospho-Smad3 (Ser425) Rabbit Polyclonal Antibody
产品货号: APRab00823



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Smad3 (Ser425) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMAD3
别名 (Alternative Names)	SMAD3; MADH3; Mothers against decapentaplegic homolog 3; MAD homolog 3; Mad3; Mothers against DPP homolog 3; hMAD-3; JV15-2; SMAD family member 3; SMAD 3; Smad3; hSMAD3
基因 ID (Gene ID)	4088
蛋白 ID (SwissProt ID)	P84022.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 58 kDa

产品名称: Phospho-Smad3 (Ser425) Rabbit Polyclonal Antibody
产品货号: APRab00823




研究背景 (Background)

Smad3 transcription factor phosphorylated and activated by TGF-beta-type receptors. A receptor-regulated Smad (R-smad). Binds directly to consensus DNA-binding elements in the promoters of target genes. In mouse required for establishment of the mucosal immune response and proper development of skeleton.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Phospho-Smad3 (Ser425) in various lysates using Phospho-Smad3 (Ser425) antibody.

注意事项 (Note)

For research use only .