

产品名称: Phospho-JAK3 (Tyr785) Rabbit Polyclonal Antibody
产品货号: APRab00826



产品概述 (Summary)

产品名称 (Production Name)	Phospho-JAK3 (Tyr785) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	JAK3
别名 (Alternative Names)	JAK3; Tyrosine-protein kinase JAK3; Janus kinase 3; JAK-3; Leukocyte janus kinase; L-JAK
基因 ID (Gene ID)	3718
蛋白 ID (SwissProt ID)	P52333.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 125 kDa; Observed MW: 125 kDa

产品名称: Phospho-JAK3 (Tyr785) Rabbit Polyclonal Antibody
产品货号: APRab00826




研究背景 (Background)

Tyrosine kinase of the non-receptor type, involved in the interleukin-2 and interleukin-4 signaling pathway. Phosphorylates STAT6, IRS1, IRS2 and PI3K.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Western blot analysis of Phospho-JAK3 (Tyr785) in 453 lysates using Phospho-JAK3 (Tyr785) antibody.




EnzymeLinked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phospho-peptide (Phospho-left) and NonPhospho-peptide (Phospho-right), using JAK3 (Phospho-Tyr785) antibody.

产品名称: Phospho-JAK3 (Tyr785) Rabbit Polyclonal Antibody
产品货号: APRab00826



Immunohistochemistry analysis of paraffin-embedded Human brain using JAK3 (Phospho-Tyr785) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Phospho-JAK3 (Tyr785) in Jurkat lysates using Phospho-JAK3 (Tyr785) antibody. The lane on the right is blocked with the Phospho- peptide.

注意事项 (Note)

For research use only .