

产品名称: Methyl-Lysine Rabbit Polyclonal Antibody  
产品货号: APRab00846



## 产品概述 (Summary)

|                        |                                          |
|------------------------|------------------------------------------|
| 产品名称 (Production Name) | Methyl-Lysine Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody               |
| 宿主 (Host)              | Rabbit                                   |
| 应用 (Application)       | WB                                       |
| 种属反应性 (Reactivity)     | Human, Mouse, Rat                        |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Methylated                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                            |
| 克隆 (Clonality)      | Polyclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Affinity Purification                                                                          |

## 免疫原信息 (Immunogen)

|                        |   |
|------------------------|---|
| 基因名 (Gene Name)        | - |
| 别名 (Alternative Names) | - |
| 基因 ID (Gene ID)        | - |
| 蛋白 ID (SwissProt ID)   | - |

## 产品应用 (Application)

|                          |                  |
|--------------------------|------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000  |
| 蛋白分子量 (Molecular Weight) | Refer to figures |

## 研究背景 (Background)

Histone lysine methylation is regarded as a very stable modification with important functions in epigenetic gene control

产品名称: Methyl-Lysine Rabbit Polyclonal Antibody  
产品货号: APRab00846




and for organizing chromatin domains.

## 研究领域 (Research Area)

Epigenetics and Nuclear Signaling

## 图片 (Image Data)



Western blot analysis of Methyl-Lysine in HeLa, 3T3, rat Brain lysates using Methyl-Lysine antibody.

## 注意事项 (Note)

For research use only .