

产品名称: Acetyl-Histone H3 (Lys9) Rabbit Polyclonal Antibody
产品货号: APRab00847



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-Histone H3 (Lys9) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H3C1
别名 (Alternative Names)	H3K9ac; H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FJ; H3C10; H3C11; HIST1H3J
基因 ID (Gene ID)	8350
蛋白 ID (SwissProt ID)	P68431.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15-17 kDa

产品名称: Acetyl-Histone H3 (Lys9) Rabbit Polyclonal Antibody
产品货号: APRab00847




研究背景 (Background)

H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Acetyl-Histone H3 (Lys9) in HeLa lysates untreated or treated with TSA lysates using Acetyl-Histone H3 (Lys9) antibody(upper) and Histone H3 (1G1) antibody(lower) TDY055.

注意事项 (Note)

For research use only .