

产品名称: Acetyl-Lysine Rabbit Polyclonal Antibody
产品货号: APRab00848



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-Lysine Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	-
别名 (Alternative Names)	-
基因 ID (Gene ID)	-
蛋白 ID (SwissProt ID)	-

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Refer to figures

研究背景 (Background)

Acetylation of lysine residues within proteins has emerged as an important mechanism used by cells to overcome this

产品名称: Acetyl-Lysine Rabbit Polyclonal Antibody
产品货号: APRab00848




repression.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Acetyl-Lysine in HeLa, 3T3, rat Brain lysates using Acetyl Lysine antibody.

注意事项 (Note)

For research use only .