

产品名称: Phospho-NMDAR2A/2B (Tyr1246/Tyr1252) Rabbit
Polyclonal Antibody
产品货号: APRab00906



产品概述 (Summary)

产品名称 (Production Name)	Phospho-NMDAR2A/2B (Tyr1246/Tyr1252) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRIN2A/GRIN2B
别名 (Alternative Names)	GRIN2A; NMDAR2A; Glutamate [NMDA] receptor subunit epsilon-1; N-methyl D-aspartate receptor subtype 2A; NMDAR2A; NR2A; hNR2A; GRIN2B; NMDAR2B; Glutamate [NMDA] receptor subunit epsilon-2; N-methyl D-aspartate receptor subtype 2B; NMDAR2B; N
基因 ID (Gene ID)	2903/2904
蛋白 ID (SwissProt ID)	Q12879/Q13224.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
----------------------	---

产品名称: Phospho-NMDAR2A/2B (Tyr1246/Tyr1252) Rabbit Polyclonal Antibody
产品货号: APRab00906

蛋白分子量 (Molecular Weight) -


研究背景 (Background)

NMDA receptor subtype of glutamate-gated ion channels possesses high calcium permeability and voltage-dependent sensitivity to magnesium. Activation requires binding of agonist to both types of subunits.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human brain using Phospho-NMDAR2A/2B (Tyr1246/Tyr1252) antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Sample with blocking peptide on the right.

注意事项 (Note)

For research use only .