

产品名称: Phospho-Annexin 2 (Ser26) Rabbit Polyclonal Antibody
产品货号: APRab00929



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Annexin 2 (Ser26) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANXA2
别名 (Alternative Names)	ANXA2; ANX2; ANX2L4; CAL1H; LPC2D; Annexin A2; Annexin II; Annexin-2; Calpactin I heavy chain; Calpactin-1 heavy chain; Chromobindin-8; Lipocortin II; Placental anticoagulant protein IV; PAP-IV; Protein I; p36
基因 ID (Gene ID)	302
蛋白 ID (SwissProt ID)	P07355.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 39 kDa

产品名称: Phospho-Annexin 2 (Ser26) Rabbit Polyclonal Antibody
产品货号: AP Rab00929




研究背景 (Background)

Originally identified as a protein inhibitor of phospholipase A2, annexin A2 has subsequently been shown to interact with an array of protein and non-protein partners, including F-actin, spectrin, SNARE complexes, RNA, and virus particles.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Phospho-Annexin 2 (Ser26) in HepG2 lysates using Phospho-ANXA2 (Ser26) antibody.

注意事项 (Note)

For research use only .