

产品名称: Phospho-XIAP (Ser87) Rabbit Polyclonal Antibody
产品货号: APRab00936



产品概述 (Summary)

产品名称 (Production Name)	Phospho-XIAP (Ser87) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	XIAP
别名 (Alternative Names)	API3; ILP1; MIHA; XLP2; BIRC4; IAP-3; hIAP3; hIAP-3
基因 ID (Gene ID)	331
蛋白 ID (SwissProt ID)	P98170.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa; Observed MW: 57 kDa

研究背景 (Background)

This gene encodes a protein that belongs to a family of apoptotic suppressor proteins. Members of this family share a

产品名称: Phospho-XIAP (Ser87) Rabbit Polyclonal Antibody
产品货号: APRab00936




conserved motif termed, baculovirus IAP repeat, which is necessary for their anti-apoptotic function. This protein functions through binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2 and inhibits apoptosis induced by menadione, a potent inducer of free radicals, and interleukin 1-beta converting enzyme.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Phospho-XIAP (Ser87) in HepG2 lysates using Phospho-XIAP (Ser87) antibody.

注意事项 (Note)

For research use only .