

产品名称: Phospho-Syk (Tyr348) Rabbit Polyclonal Antibody
产品货号: APRab00940



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Syk (Tyr348) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SYK
别名 (Alternative Names)	SYK; Tyrosine-protein kinase SYK; Spleen tyrosine kinase; p72-Syk
基因 ID (Gene ID)	6850
蛋白 ID (SwissProt ID)	P43405.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 72 kDa; Observed MW: 72 kDa

研究背景 (Background)

This gene encodes a member of the family of non-receptor type Tyr protein kinases. This protein is widely expressed in


产品名称: Phospho-Syk (Tyr348) Rabbit Polyclonal Antibody
产品货号: APRab00940

hematopoietic cells and is involved in coupling activated immunoreceptors to downstream signaling events that mediate diverse cellular responses, including proliferation, differentiation, and phagocytosis. It is thought to be a modulator of epithelial cell growth and a potential tumour suppressor in human breast carcinomas. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Phospho-Syk (Tyr348) in 3T3 lysates using Phospho-Syk (Tyr348) antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Phospho-Syk (Tyr348) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.

注意事项 (Note)

For research use only .