

产品名称: Phospho-PYK2 (Tyr402) Rabbit Polyclonal Antibody
产品货号: APRab00952



产品概述 (Summary)

产品名称 (Production Name)	Phospho-PYK2 (Tyr402) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PTK2B
别名 (Alternative Names)	PTK2B; FAK2; PYK2; RAFTK; Protein-tyrosine kinase 2-beta; Calcium-dependent tyrosine kinase; CADTK; Calcium-regulated non-receptor proline-rich tyrosine kinase; Cell adhesion kinase beta; CAK-beta; CAKB; Focal adhesion kinase 2; FADK 2; Pro
基因 ID (Gene ID)	2185
蛋白 ID (SwissProt ID)	Q14289.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:100-1:200
----------------------	---

产品名称: Phospho-PYK2 (Tyr402) Rabbit Polyclonal Antibody
产品货号: APRab00952



蛋白分子量 (Molecular Weight) Calculated MW: 116 kDa; Observed MW: 116 kDa


研究背景 (Background)

PYK2 a nonreceptor tyrosine kinase of the Fak family. Predominantly expressed in the cells derived from hematopoietic lineages and in the central nervous system. PYK2 is one of the signaling mediators for G-protein-coupled receptors. Involved in calcium induced regulation of ion channel and activation of the map kinase signaling pathway. Interacts with the SH2 domain of Grb2.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Western blot analysis of Phospho-PYK2 (Tyr402) in 293 lysates using Phospho-PYK2 (Tyr402) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using Pyk2(Phospho-Tyr402) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.

产品名称: Phospho-PYK2 (Tyr402) Rabbit Polyclonal Antibody
产品货号: APRab00952



Immunofluorescence analysis of Phospho-PYK2 (Tyr402) in Hela cells using Pyk2(Phospho-Tyr402) antibody(red).

注意事项 (Note)

For research use only .