

产品名称: Phospho-YAP (Ser127) Rabbit Polyclonal Antibody
产品货号: APRab00954



产品概述 (Summary)

产品名称 (Production Name)	Phospho-YAP (Ser127) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YAP1
别名 (Alternative Names)	YAP1; YAP65; Yorkie homolog; 65 kDa Yes-associated protein; YAP65
基因 ID (Gene ID)	10413
蛋白 ID (SwissProt ID)	P46937.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 65 kDa

研究背景 (Background)

YAP (Yes-associated protein, YAP65) was identified based on its ability to associate with the SH3 domain of Yes.

产品名称: Phospho-YAP (Ser127) Rabbit Polyclonal Antibody
产品货号: APRab00954




Transcriptional regulator which can act both as a coactivator and a corepressor and is the critical downstream regulatory target in the Hippo signaling pathway that plays a pivotal role in organ size control and tumor suppression by restricting proliferation and promoting apoptosis.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of Phospho-YAP (Ser127) in COS7 lysates using Phospho-YAP (Ser127) antibody.



Immunohistochemistry analysis of paraffin-embedded Human heart using YAP (Phospho-Ser127) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



产品名称: Phospho-YAP (Ser127) Rabbit Polyclonal Antibody
产品货号: APRab00954



Immunofluorescence analysis of Phospho-YAP (Ser127) in Hela cells using YAP (Phospho- Ser127) antibody(green) and beta-tubulin antibody(green).

注意事项 (Note)

For research use only .