

产品名称: Phospho-ALK (Tyr1604) Rabbit Polyclonal Antibody
产品货号: APRab00955



产品概述 (Summary)

产品名称 (Production Name)	Phospho-ALK (Tyr1604) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALK
别名 (Alternative Names)	ALK; ALK tyrosine kinase receptor; Anaplastic lymphoma kinase; CD antigen CD246
基因 ID (Gene ID)	238
蛋白 ID (SwissProt ID)	Q9UM73.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 176 kDa; Observed MW: 176 kDa

产品名称: Phospho-ALK (Tyr1604) Rabbit Polyclonal Antibody
产品货号: APRab00955




研究背景 (Background)

Involved in diverse cellular processes such as ribosome biogenesis, centrosome duplication, protein chaperoning, histone assembly, cell proliferation, and regulation of tumor suppressors p53/TP53 and ARF. Binds ribosome presumably to drive ribosome nuclear export. Associated with nucleolar ribonucleoprotein structures and bind single-stranded nucleic acids.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)




Western blot analysis of Phospho-ALK (Tyr1604) in HepG2 lysates using Phospho-ALK (Tyr1604) antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain tissue, using ALK (Phospho-Tyr1604) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.

产品名称: Phospho-ALK (Tyr1604) Rabbit Polyclonal Antibody
产品货号: APRab00955



Immunofluorescence analysis of Phospho-ALK (Tyr1604) in HeLa cells using ALK (Phospho-Tyr1604) antibody(green).

注意事项 (Note)

For research use only .