

产品名称: Annexin V Rabbit Polyclonal Antibody (HRP)
产品货号: APRab01310



产品概述 (Summary)

产品名称 (Production Name)	Annexin V Rabbit Polyclonal Antibody (HRP)
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	HRP
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANXA5
别名 (Alternative Names)	PP4; ANX5; ENX2; RPRGL3.
基因 ID (Gene ID)	308
蛋白 ID (SwissProt ID)	P08758.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 36 kDa; Observed MW: 30 kDa

研究背景 (Background)

This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is

产品名称: Annexin V Rabbit Polyclonal Antibody (HRP)
产品货号: AP Rab01310




involved in the blood coagulation cascade.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Annexin V in HeLa and COS7 lysates using Annexin V antibody.conjugated with HRP

注意事项 (Note)

For research use only .