

产品名称: Phospho-PKA RII alpha (Ser99) Rabbit Polyclonal Antibody
产品货号: APRab01319



产品概述 (Summary)

产品名称 (Production Name)	Phospho-PKA RII alpha (Ser99) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Other

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRKAR2A
别名 (Alternative Names)	PRKAR2A; PKR2; PRKAR2; cAMP-dependent protein kinase type II-alpha regulatory subunit
基因 ID (Gene ID)	5576
蛋白 ID (SwissProt ID)	P13861.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 46 kDa; Observed MW: 51 kDa

产品名称: Phospho-PKA RII alpha (Ser99) Rabbit Polyclonal Antibody
产品货号: APRab01319




研究背景 (Background)

Regulatory subunit of the cAMP-dependent protein kinases involved in cAMP signaling in cells. Type II regulatory chains mediate membrane association by binding to anchoring proteins, including the MAP2 kinase.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Phospho-PKA R2 (Ser99) in K562 lysates using Phospho-PKA RII alpha (Ser99) antibody.

注意事项 (Note)

For research use only .