

产品名称: FUBP1 Rabbit Polyclonal Antibody
产品货号: APRab01321



产品概述 (Summary)

产品名称 (Production Name)	FUBP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	FUBP1
别名 (Alternative Names)	DNA helicase V; FBP; FUBP; Fubp1; Fubp4; FUSE-binding protein 1; HDH V
基因 ID (Gene ID)	8880
蛋白 ID (SwissProt ID)	Q96AE4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 74 kDa

研究背景 (Background)

Regulates MYC expression by binding to a single-stranded far-upstream element upstream of the MYC promoter. May act

产品名称: FUBP1 Rabbit Polyclonal Antibody
产品货号: APRab01321




both as activator and repressor of transcription.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of FUBP1 in Jurkat, K562 and HeLa lysates using FUBP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human spleen using FUBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .