

产品名称: Acetyl-p53 (Lys370) Rabbit Polyclonal Antibody
产品货号: APRab01331



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-p53 (Lys370) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	TP53
别名 (Alternative Names)	Antigen NY-CO-13; Cellular tumor antigen p53; Phosphoprotein p53; TP53; Tumor suppressor p53
基因 ID (Gene ID)	7157
蛋白 ID (SwissProt ID)	P04637.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 53 kDa

研究背景 (Background)

产品名称: Acetyl-p53 (Lys370) Rabbit Polyclonal Antibody
产品货号: APRab01331




Tumor protein p53, a nuclear protein, plays an essential role in the regulation of cell cycle, specifically in the transition from G0 to G1. It is found in very low levels in normal cells, however, in a variety of transformed cell lines, it is expressed in high amounts, and believed to contribute to transformation and malignancy. p53 is a DNA-binding protein containing DNA-binding, oligomerization and transcription activation domains.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of p53 (acetyl K370) in HeLa lysates treated with Trichostatin A using Acetyl-p53 (Lys370) antibody.

注意事项 (Note)

For research use only .