

产品名称: CD8 alpha Rabbit Polyclonal Antibody
产品货号: APRab01333



产品概述 (Summary)

产品名称 (Production Name)	CD8 alpha Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD8A
别名 (Alternative Names)	CD8A; MAL; T-cell surface glycoprotein CD8 alpha chain; T-lymphocyte differentiation antigen T8/Leu-2; CD antigen CD8a
基因 ID (Gene ID)	925
蛋白 ID (SwissProt ID)	P01732.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 32 kDa

研究背景 (Background)

产品名称: CD8 alpha Rabbit Polyclonal Antibody
产品货号: APRab01333




Identifies cytotoxic/suppressor T-cells that interact with MHC class I bearing targets. CD8 is thought to play a role in the process of T-cell mediated killing. CD8 alpha chains binds to class I MHC molecules alpha-3 domains.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of CD8 in Jurkat lysates using CD8 alpha antibody.




Immunohistochemistry analysis of paraffin-embedded Human colon using CD8 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of CD8 alpha in Jurkat using CD8 antibody.

产品名称: CD8 alpha Rabbit Polyclonal Antibody
产品货号: APRab01333



Western blot analysis of CD8 alpha in Raji lysates using CD8 antibody.

注意事项 (Note)

For research use only .