

产品名称: PKC delta Rabbit Polyclonal Antibody
产品货号: APRab01337



产品概述 (Summary)

产品名称 (Production Name)	PKC delta Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRKCD
别名 (Alternative Names)	Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta
基因 ID (Gene ID)	5580
蛋白 ID (SwissProt ID)	Q05655.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 78 kDa; Observed MW: 78 kDa

研究背景 (Background)

Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the

产品名称: PKC delta Rabbit Polyclonal Antibody
产品货号: AP Rab01337




second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of PKC delta in HeLa lysates using PKC delta antibody.

注意事项 (Note)

For research use only .