

产品名称: Human IgM Rabbit Polyclonal Antibody
产品货号: APRab01340



产品概述 (Summary)

产品名称 (Production Name)	Human IgM Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGHM
别名 (Alternative Names)	Ig mu chain C region
基因 ID (Gene ID)	3507.0
蛋白 ID (SwissProt ID)	P01871.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 49 kDa; Observed MW: 75 kDa

研究背景 (Background)

IgM measurement yields information about the body's immediate resistance and response to infection as well as

产品名称: Human IgM Rabbit Polyclonal Antibody
产品货号: APRab01340




information related to specific diseases. Decreased levels are associated with immune deficiency states, hereditary deficiencies, and myeloma. Increased levels can be associated with Waldenstrom's macroglobulinemia, chronic infection and hepatocellular disease.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of Human IgM in Human plasma lysates using Human IgM antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Human IgM antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .