

**产品名称:** Angiotensin II Type 1 Receptor Rabbit Polyclonal Antibody  
**产品货号:** APRab01341



## 产品概述 (Summary)

|                        |                                                           |
|------------------------|-----------------------------------------------------------|
| 产品名称 (Production Name) | Angiotensin II Type 1 Receptor Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody                                |
| 宿主 (Host)              | Rabbit                                                    |
| 应用 (Application)       | WB                                                        |
| 种属反应性 (Reactivity)     | Human, Mouse, Rat                                         |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Chromatography                                                                           |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                         |
|------------------------|-------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | AGTR1                                                                                                                   |
| 别名 (Alternative Names) | AGTR1; AGTR1A; AGTR1B; AT2R1; AT2R1B; Type-1 angiotensin II receptor; AT1AR; AT1BR; Angiotensin II type-1 receptor; AT1 |
| 基因 ID (Gene ID)        | 185                                                                                                                     |
| 蛋白 ID (SwissProt ID)   | P30556.                                                                                                                 |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000                            |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 41 kDa; Observed MW: 41 kDa |

产品名称: Angiotensin II Type 1 Receptor Rabbit Polyclonal  
Antibody  
产品货号: AP Rab01341




## 研究背景 (Background)

Receptor for angiotensin II. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

## 研究领域 (Research Area)

Cancer

## 图片 (Image Data)



Western blot analysis of AGTR1 in (1) HeLa lysates; (2) PC-12 lysates using Angiotensin II Type 1 Receptor antibody.

## 注意事项 (Note)

For research use only .