

产品名称: MSR1 Rabbit Polyclonal Antibody
产品货号: APRab01349



产品概述 (Summary)

产品名称 (Production Name)	MSR1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	MSR1
别名 (Alternative Names)	Macrophage scavenger receptor types I and II; Macrophage acetylated LDL receptor I and II; Scavenger receptor class A member 1; CD204; MSR1; SCARA1
基因 ID (Gene ID)	4481
蛋白 ID (SwissProt ID)	P21757.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 75 kDa

研究背景 (Background)

产品名称: MSR1 Rabbit Polyclonal Antibody
产品货号: APRab01349




Membrane glycoproteins implicated in the pathologic deposition of cholesterol in arterial walls during atherogenesis. Two types of receptor subunits exist. These receptors mediate the endocytosis of a diverse group of macromolecules, including modified low density lipoproteins (LDL). Isoform III does not internalize actetylated LDL.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of MSR1 in (1) HeLa lysates; (2) mouse placenta lysates using MSR1 antibody.

注意事项 (Note)

For research use only .