

产品名称: RBBP4 Rabbit Polyclonal Antibody
产品货号: APRab01353



产品概述 (Summary)

产品名称 (Production Name)	RBBP4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBBP4
别名 (Alternative Names)	NURF55; RBAP48; lin-53
基因 ID (Gene ID)	5928
蛋白 ID (SwissProt ID)	Q09028.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 48 kDa

研究背景 (Background)

Core histone-binding subunit that may target chromatin assembly factors, chromatin remodeling factors and histone

产品名称: RBBP4 Rabbit Polyclonal Antibody
产品货号: APRab01353




deacetylases to their histone substrates in a manner that is regulated by nucleosomal DNA. Component of several complexes which regulate chromatin metabolism.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of RbAp48 in Jurkat lysates using RBBP4 antibody.



Western blot analysis of RBBP4 in 3T3, HCT116, COS7 and HeLa lysates using RbAp48 antibody.

注意事项 (Note)

For research use only .