

产品名称: ATP6V1E1 Rabbit Polyclonal Antibody
产品货号: APRab01359



产品概述 (Summary)

产品名称 (Production Name)	ATP6V1E1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATP6V1E1
别名 (Alternative Names)	V-ATPase subunit E 1; p31
基因 ID (Gene ID)	529
蛋白 ID (SwissProt ID)	P36543.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 26 kDa

研究背景 (Background)

Subunit of the V1 complex of vacuolar(H⁺)-ATPase (V-ATPase), a multisubunit enzyme composed of a peripheral complex

产品名称: ATP6V1E1 Rabbit Polyclonal Antibody
产品货号: APRab01359




(V1) that hydrolyzes ATP and a membrane integral complex (V0) that translocates protons.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of ATP6V1E1 in 3T3, HeLa, 293 lysates using ATP6V1E1 antibody.



Immunocytochemistry analysis of ATP6V1E1(green) in HeLa using ATP6V1E1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ATP6V1E1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: ATP6V1E1 Rabbit Polyclonal Antibody
产品货号: APRab01359



注意事项 (Note)

For research use only .