

产品名称: GDF15 Rabbit Polyclonal Antibody
产品货号: APRab01365



产品概述 (Summary)

产品名称 (Production Name)	GDF15 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GDF15 GDF15; MIC1; PDF; PLAB; PTGFB; Growth/differentiation factor 15; GDF-15;
别名 (Alternative Names)	Macrophage inhibitory cytokine 1; MIC-1; NSAID-activated gene 1 protein; NAG-1; NSAID-regulated gene 1 protein; NRG-1; Placental TGF-beta; Placental bone morphogenetic
基因 ID (Gene ID)	9518
蛋白 ID (SwissProt ID)	Q99988.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 35,13 kDa

产品名称: GDF15 Rabbit Polyclonal Antibody
产品货号: APRab01365




研究背景 (Background)

Highly expressed in placenta, with lower levels in prostate and colon and some expression in kidney.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of GDF15/Mic1 in HeLa, MCF-7 lysates using GDF15 antibody.

注意事项 (Note)

For research use only .